## Spezialpunktkanüle DPK 2035

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## Special puncture needle DPK 2035

### COMPONENTS

Special non-coring needle for port puncture with cannula guide pin and needle-free adapter (NFA) for implanted port catheter systems

External diameter of 2.0 mm, internal diameter 1.80 mm, internal volume 0.27 ml, gauge: 14 G

### INDICATIONS

Access to implanted portcatheter systems (e.g. 111252 D, 111255 D, 111263 APH, 111259 AS/PT) e.g. for hemodialysis, peritoneal dialysis, peritoneal chemotherapy, extracorporeal apheresis, ascites therapy, and drainage of fluids. Use only the appropriate non-coring special puncture needles to puncture the silicone membranes of the port chambers. These are specially beveled to prevent silicone particle fragmentation and membrane puncture defects.

In order to ensure an adequate volume (e.g. of blood for hemodialysis), large bore cannulas are required (e.g. DPK 2035).

Flow approx 500 ml/min (hemodialysis), gravitation 300 ml/min

**Patient group:** Patients with the need for long-term central access to the vascular system for the above applications.

**Operator:** Medical professionals (doctors, nurses) with appropriate training.

### CONTRAINDICATIONS

- Inflammation or infection at portal site
- Coagulation disorders
- Portal system leakage

### WARNINGS

Due to incorrect use during puncture of the port air can be sucked into the system via the needle. Therefore the closure mechanism of the needle/infusion tubing should always be closed during puncture or when changing syringes/infusions.

**Cave:** Puncture of the port should be performed in the supine position to prevent air embolism (in case of vascular access). Negative pressure gradients develop in the superior vena cava in the sitting position or when the torso is elevated. Patient's upper body is upright.

### CARE OF THE SYSTEM / SPECIAL APPLICATION NOTES

On each occasion before treatment, ensure the correct position of the port chamber by palpation, and check that there is no wound, hematoma or infection present.

General standards in hygiene for handling of portcatheter systems must be observed.

Removal of the cannula after use is recommended.

**Cave:** Blood deposits in the system generally indicate faulty handling, leakage or a defective septum inside the port

### INSTRUCTIONS FOR USE

Check proper position of cannula guide pin before use.

- have patient in supine position (if possible)
- skin disinfection, sterile gloves, mask (patient and operator)
- palpation of silicone membrane of the port
- insert special puncture needle (with cannula guide pin) through skin and tissues through the silicone membrane until tip of needle touches port wall.

**Cave:** air embolism

- remove cannula guide pin (NFA keeps the system closed)
- if necessary attach a syringe to the needle-free adapter and aspirate residual volume content from port chamber and catheter (approx. 3 ml)
- check patency of system with e.g. physiologic saline/dialyzing fluid and attach therapy system
- upon termination of therapy or after and between injections of pharmaceuticals flush and block the port chamber if necessary (flush

with saline solution 0,9%, block with e.g. heparin, citrate or taurolidine as instructed by the physician)

- remove the puncture needle under slight pressure, in order to prevent intrusion of body fluid into the system or development of underflow
- skin disinfection, cover dry

Use only the appropriate non-coring special puncture cannula for puncturing - use the entire membrane (see fig. back cover)

Do not use syringes < 10 ml (danger of overload pressure)

Catheter or port obstruction due to thrombosis can generally be rectified by lysis with streptokinase/ urokinase (e.g. 5000 IU/ml in 0.9% NaCl, by injecting small amounts of the solution at regular intervals; allow the procedure to take effect and check patency). Previous: x-ray control

Potential reasons for insufficient or failing flow:

- The tip of the needle is not positioned correctly in the port cavity
- Fibrin sedimentation or a thrombus occluded catheter or port
- The catheter is kinked
- The tip of the catheter is not placed correctly
- If blood aspiration is not successful, the port chamber should be initially rinsed and the patient's position changed slightly.

The concentrations given are recommendations and must always be adjusted to the general situation of each patient. The treating doctor is responsible for this. For all the drugs and liquids listed, we refer primarily to the instruction leaflets provided by the respective pharmaceutical manufacturers. The specific features of therapy have to be taken into account. Pay attention to the instructions for use provided with the port system. The port systems are merely access systems with no effect on function.

### PRECAUTIONS

Meticulous hygiene and sterile handling technique are mandatory for this system, e.g. Puncture to be carried out only under sterile conditions

Disinfect hands and skin before using the system. Use sterile products only

The medical product may only be used once; it is not suitable for reprocessing. In case of reuse, complete functionality cannot be guaranteed.

### MONITORING

Continuous monitoring is ensured per se, since the system must be used

and subsequently maintained on a regular basis. The system should not be used by doctors/medical staff who are unfamiliar with the product and/or potential complications. Complications may occur at any time, both during and after the invasive procedure.


### DISPOSAL

Please use needle drop containers for disposal to avoid injury from needle sticks or contamination from body fluids.

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Special puncture needle  
 Ref.: DPK 2035

Instructions for Use  
 please read carefully!

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