Spezialpunktkanüle SFN 0930 B

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SFN - G port needle (straight)

COMPONENTS

Special non-coring needle with luer-lock connector for implanted port catheter systems

SFN 0930 B: outer diameter 0.9 mm, inner diameter 0.6 mm, length: 30 mm, inner volume approx. 10 ml (with tubing), gauge: 20 G

The **connecting tube** is made of polyurethane (PU): Ø 1.6 mm; length: 130 mm.

INDICATIONS

Access to implanted portcatheter systems e. g. for rinsing of the system, administration of pharmaceuticals or infusions or blood sampling.

Use only the appropriate non-coring special puncture needles to puncture the silicone membranes of the port chambers. These are specially beveled to prevent silicone particle fragmentation and membrane puncture defects.

A Luer-Lock connector allows connection to syringes, infusion systems, etc. with a 6% Luer taper for intravascular and hypodermic applications.

Patient group: Patients with the need for long-term central access to the vascular system for the above applications.

Operator: Medical professionals (doctors, nurses) with appropriate training.

CONTRAINDICATIONS

- Inflammation or infection at portal site
- Coagulation disorders
- Portal system leakage

WARNINGS

Due to incorrect use during puncture of the port air can be sucked into the system via the needle. Therefore the closure mechanism of the needle/infusion tubing should always be closed during puncture or when changing syringes/infusions.

Cave: Puncture of the port should be performed in the supine position to prevent air embolism. Negative pressure gradients develop in the superior vena cava in the sitting position or when the torso is elevated.

CARE OF THE SYSTEM / SPECIAL APPLICATION NOTES

On each occasion before treatment, ensure the correct position of the port chamber by palpation, and check that there is no wound or infection present.

General standards in hygiene for handling of portcatheter systems must be observed. Removal of the cannula after use is recommended.

Cave: Blood deposits in the system generally indicate faulty handling, leakage or a defective septum inside the port.

INSTRUCTIONS FOR USE

- have patient in supine position (if possible)
- skin disinfection, sterile gloves, mask (patient and operator)
- remove air from needle e.g. by flushing with saline
- please note that needle and tubing (a procedure pack) are pre-connected via Luer-Lock and the connection should be reconfirmed before use (tightness and patency).
- palpation of silicone membrane of the port
- advance needle quickly through skin and tissues through the silicone membrane until tip of needle touches the metallic base of the port
- if necessary aspirate residual volume content from portchamber and catheter (approx. 3 ml)
- check patency of system with e.g. physiologic saline/dialyzing fluid and attach therapy system
- upon termination of therapy or after and between injections of pharmaceuticals flush and block the port chamber if necessary (flush with saline solution 0.9 %, block with e. g. heparin, citrate or taurolidin as instructed by the physician)
- remove the puncture needle under slight pressure, in order to prevent intrusion of blood into the system or development of underflow
- skin disinfection, cover dry

No danger of bending of the needle tip due to special bevel of the tip. Use only the appropriate non-coring special puncture cannula for puncturing - use the entire membrane.

Do not use syringes < 10 ml (danger of overload pressure).

Catheter or port obstruction due to thrombosis can generally be rectified by lysing with streptokinase/urokinase (e. g. 5000 IU/ml in 0.9 % NaCl, by injecting small amounts of the solution at regular intervals; allow the procedure to take effect and check patency). Previous: x-ray control.

Potential reasons for insufficient or failing flow:

- The tip of the needle is not positioned correctly in the port cavity
- Fibrin sedimentation or a thrombus occluded catheter or port
- The catheter is kinked
- The tip of the catheter is not placed correctly
- If blood aspiration is not successful, the port chamber should be initially rinsed and the patient's position changed slightly.

The concentrations given are recommendations and must always be adjusted to the general situation of each patient. The treating doctor is responsible for this. For all the drugs and liquids listed, we refer primarily to the instruction leaflets provided by the respective pharmaceutical manufacturers. The specific features of therapy have to be taken into account. Pay attention to the instructions for use provided with the port system. The port systems are merely access systems with no effect on function.

PRECAUTIONS

Meticulous hygiene and sterile handling technique are mandatory for this system, e.g.

- Puncture to be carried out only under sterile conditions
- Disinfect hands and skin before using the system
- Use sterile products only! If the sterile packaging was already opened or damaged, DO NOT use the product!

The medical product may only be used once; it is not suitable for reprocessing. In case of reuse, complete functionality cannot be guaranteed.

MONITORING

Continuous monitoring is ensured per se, since the system must be used and subsequently maintained on a regular basis. The system should not be used by doctors/medical staff who are unfamiliar with the product and/or potential complications. Complications may occur at any time, both during and after the invasive procedure.

The medical device is for single use only. It is not appropriate for reprocessing.

DISPOSAL

Please use needle drop containers for disposal to avoid injury from needle sticks or contamination from body fluids.

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Notice to operators and patients: Please report all serious incidents related to this product to the manufacturer, as well as to the competent authority of the Member State where you are located.

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
Portneedle
Ref.: SFN 0930 B - straight needle
with tube

Instructions for Use
please read carefully!

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