## Spezialpunktkanüle SFN 1835 B, SFN 1535 B

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## Special puncture needle SFN 1835 B, SFN 1535 B

### COMPONENTS

Special non-coring needle with luer-lock connector grip and connective tube with clamp for implanted port catheter systems

**SFN 1835 B:** External diameter of 1.80 mm, internal diameter of 1.50 mm, internal volume 10 ml (with tubing), gauge: 15 G

**SFN 1535 B:** External diameter of 1.50 mm, internal diameter of 1.20 mm, internal volume 10 ml (with tubing), gauge: 17 G

The **connecting tube** is made of DEHP-free polyvinyl chloride (PVC): Ø 3.5 mm; length: 100 mm.

### INDICATIONS

Access to implanted portcatheter systems (e.g. 111252 D, 111255 D, 111259 AS/PT, 111262 DR, 111263 APH, e.g. for hemodialysis, peritoneal dialysis, peritoneal chemotherapy, apheresis or ascites therapy).

Use only the appropriate non-coring special puncture needles to puncture the silicone membranes of the port chambers. These are specially beveled to prevent silicone particle fragmentation and membrane puncture defects.

A Luer-Lock connector allows connection to syringes, infusion systems, etc. with a 6% Luer taper for intravascular and hypodermic applications.

Flow rate of the needle approx.: (depending on the system used)

SFN 1835 B - 300 ml/min (hemodialysis), 225 ml/min (gravitation)

SFN 1535 B - 200 ml/min (hemodialysis), 150 ml/min (gravitation)

**Patient group:** Patients with the need for long-term central access to the vascular system for the above applications.

**Operator:** Medical professionals (doctors, nurses) with appropriate training.

### CONTRAINDICATIONS

- Inflammation or infection at portal site
- Coagulation disorders
- Portal system leakage

### WARNINGS

Due to incorrect use during puncture of the port air can be sucked into the system via the needle. Therefore the dosing mechanism of the needle / infusion tubing should always be dosed during puncture or when changing syringes/infusions.

**Cave:** Puncture of the port should be performed in the supine position to prevent air embolism. Negative pressure gradients develop in the superior vena cava in the sitting position or when the torso is elevated.

### CARE OF THE SYSTEM / SPECIAL APPLICATION NOTES

On each occasion before treatment, ensure the correct position of the port chamber by palpation, and check that there is no wound or infection present. General standards in hygiene for handling of portcatheter systems must be observed. Removal of the cannula after use is recommended.

**Cave:** Blood deposits in the system generally indicate faulty handling, leakage or a defective septum inside the port.

### INSTRUCTIONS FOR USE

- have patient in supine position (if possible)
- skin disinfection, sterile gloves, mask (patient and operator)
- remove air from needle and tubing by flushing with saline
- please note that needle and tubing (a procedure pack) are pre-connected via Luer-Lock and the connection should be reconfirmed before use (tightness and patency)
- palpation of silicone membrane of the port
- advance needle quickly through skin and tissues through the silicone membrane until tip of needle touches the metallic base of the port
- if necessary aspirate residual volume content from portchamber and catheter (approx. 3 ml)
- check patency of system with e.g. physiologic saline / dialyzing fluid and attach therapy system
- upon termination of therapy or after and between injections of pharmaceuticals flush and block the port chamber if necessary (flush with

saline solution 0.9 % block with e. g. heparin, citrate or taurolidine as instructed by the physician)

- remove the puncture needle under slight pressure, in order to prevent intrusion of blood into the system or development of underflow
- skin disinfection, cover dry

No danger of bending of the needle tip due to special bevel of the tip  
Use only the appropriate non-coring special puncture cannula for puncturing  
- use the entire membrane (see fig. page 4). Do not use syringes < 10 ml (danger of overload pressure).

Catheter or port obstruction due to thrombosis can generally be rectified by lysing with streptokinase/ urokinase (e. g. 5000 IU/ml in 0.9 % NaCl, by injecting small amounts of the solution at regular intervals; allow the procedure to take effect and check patency). Previous: x-ray control  
Potential reasons for insufficient or failing flow:

- The tip of the needle is not positioned correctly in the port cavity
- Fibrin sedimentation or a thrombus occluded catheter or port
- The catheter is kinked
- The tip of the catheter is not placed correctly
- If blood aspiration is not successful, the port chamber should be initially rinsed and the patient's position changed slightly.

The concentrations given are recommendations and must always be adjusted to the general situation of each patient. The treating doctor is responsible for this. For all the drugs and liquids listed, we refer primarily to the instruction leaflets provided by the respective pharmaceutical manufacturers. The specific features of therapy have to be taken into account. Pay attention to the instructions for use provided with the port system. The port systems are merely access systems with no effect on function.

### PRECAUTIONS

Meticulous hygiene and sterile handling technique are mandatory for this system, e.g.

- Puncture to be carried out only under sterile conditions
- Disinfect hands and skin before using the system
- Use sterile products only! If the sterile packaging was already opened or damaged, DO NOT use the product!

The medical product may only be used once; it is not suitable for reprocessing. In case of reuse, complete functionality cannot be guaranteed.

### MONITORING

Continuous monitoring is ensured per se, since the system must be used and subsequently maintained on a regular basis. The system should not be used by doctors/medical staff who are unfamiliar with the product and/or potential complications. Complications may occur at any time, both during and after the invasive procedure. The medical device is for single use only. It is not appropriate for reprocessing.

### DISPOSAL

Please use needle drop containers for disposal to avoid injury from needle sticks or contamination from body fluids.

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

Special puncture needle  
 Ref.: SFN 1835 B | SFN 1535 B

Instructions for Use  
 please read carefully!

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